

ambre



Ambre 8200
Mammography Chair

DESCRIPTION

A mammography exam can be an emotional experience for the patient. That is why it is important to use the right equipment that will make the exam as effective and as pleasant as possible. Some patients need to stay seated through the entire consultation. It is also easier to get the images right if the patient is seated in a stable, ergonomically correct, seated position. For more time consuming examinations the patient should be seated comfortably and relaxed.

The Ambre 8200 Mammography is an examination chair that offers the patient a correct posture and at the same time facilitates a good work environment for the staff. It also makes the examinations quicker and more cost effective.

FEATURES



- High backrest
- Neck support
- Electric forward/backward tilt via remote control
- Adjustable armrests
- White frame
- Brake lever
- Foldable footrest
- 10cm Castors



Base with footbrake and 1/4 foldable foot ring.



Push handle and wired remote control.

POSITIONS



Neutral position



Electric seat tilt



Electric backrest tilt

Ambre 8200 Mammography Chair	
Seat	44cm x 43cm
Electric seat tilt	15° forward / 15° backwards
Extra high backrest	42cm x 43cm
Armrest	25cm (alt. 30cm)
High seat height	47cm – 75cm
Base	White 48cm x 53cm
Castors	10cm XL
Brake lever	35cm
Foot ring	¼ foldable foot ring

Contact Us:

Phone: 01522 420091

Email: enquiries@ambrestone.co.uk

Website: www.ambrestone.co.uk



ENRICHING LIVES THROUGH MAXIMUM POSTURAL
SUPPORT AND PRESSURE MANAGEMENT

Ambre-stone Ltd
Fosters Yard
Station Road
Langworth
Lincoln
LN3 5BD

01522 420091
enquiries@ambrestone.co.uk
www.ambrestone.co.uk