

Thinking

EFFECTO PLUS  
**inserto**  
BACK

BACK

2 A Teschio e denti, 1489.

R.L. 19058 v (B. 41 r; K/P. 42 v).

Il vacuo della cassa dell'occhio, e 'l vacuo del-  
fosso sostenitore della guancia, e quello del  
naso e della bocca sono d'eguale profondità e  
terminano sotto il senso comune (siede centra-  
le delle sensazioni) per linea perpendicolare. E  
ciascuna d'esse vacuità ha tanto di profondità  
quanti è la terza parte del volto dell'omo, cioè  
dal mento a capegli

A NEW SOLUTION  
The Trunk Total Support



## Why Effecto Plus

 HELPS TRUNK ALIGNMENT



**Crash tested  
ISO 7176-19**



**5 modular  
configurable inserts**



**Multi-adjustable**

### CLINICAL INFORMATION

- Moderate control of the trunk
- Lateral Pelvis-Trunk instability
- Severe deformities

### TECHNICAL FEATURES

- 5 modular configurable inserts as standard
- Two layers of foam. Self-modeling polymeric foam and reticulated foam for ventilation
- Microclimatic ventilated cover
- Aluminum flat shell T shaped
- Mounting hardware kit with 4 attaching points

### SPARE PARTS

- Cover
- Mounting hardware kit

### ACCESSORIES

- Fixatis Range Harnesses
- Capitis Range Headrest
- Headrest Mounting Plate

**TTS** EMBRACING SHAPES






**inserto EFFECTO PLUS**

5 modular configurable inserts

MODEL	Paediatric			Adult				
	3033	3338	3642	3945	4245	4250	4550	4855
Width (cm)	30	33	36	39	42	42	45	48
Height (cm)	33	38	42	45	45	50	50	55

• Maximum Load (referred to size 4855): **136 KG**

• Compatible with wheelchair tubing diameter from 19 up to 25 mm. and from ¾" up to 1".



## Technical Information

Backrest Technical Information				
SIZE (cm)	Shell Height (cm)	Mounting on wheelchairs (outside tubes) (cm)	Reclination (°)	Rotation (°)
30x33	33	from 27 up to 33	-7 up to +7	-5 up to +5
33x38	38	from 30 up to 36		
36x42	42	from 33 up to 39		
39x45	45	from 36 up to 42		
42x45	45	from 39 up to 45		
42x50	50	from 39 up to 45		
45x50	50	from 42 up to 48		
48x55	55	from 45 up to 51		

### COMPATIBILITY AND REMOVAL

*It is compatible with the most common wheelchairs tubing size from 19 up to 25 mm and from ¾" up to 1".  
Quick removal.*



## Accessories



### Headrest *Versa* Capitis<sup>®</sup> HEAD SUPPORTS

#### MOUNTING SYSTEM



#### CAPITIS CONFORT



#### CAPITIS ANATOMICO




CAPITIS CONFORT DIMENSIONS (expressed in centimeters)

	A	B	C	D	E	F
	21	14	13	15	08	04
	30	14	13	24	10	05
	42	16	14	34	14	07

CAPITIS ANATOMICO DIMENSIONS (expressed in centimeters)

	A	B	C	D	E	F
	18	13	08	13	04	02
	24	17	12	19	06	04
	28	19	13	21	08	06
	32	25	16	23	08	08


### Harnesses *Versa* Fixatis<sup>®</sup> FASTENERS



ARTICLE	A	B	C	D	E
PMD 16XS	11	18	17	05	29
PMD 16S	17	20	21	05	34
PMD 16M	20	20	26	05	46
PMD 16L	20	25	30	06	52



ARTICLE	L	H
PMD 17S	14	16
PMD 17M	16	18
PMD 17L	23	23



- H = Height of the pad
- L = Width of the pad


ARTICLE	H	L
PMD 21XS	25	3.8
PMD 21S	35	3.8
PMD 21M	38	5
PMD 21L	45	5



MANUFACTURED BY:



Via Montagna Z.I. . Lotto 41  
72023 Mesagne (Br) - ITALY  
Tel . + (39) 0831. 777 840  
Fax . + (39) 0831. 730 739  
Web site: [www.promedicare.it](http://www.promedicare.it)  
e mail: [info@promedicare.it](mailto:info@promedicare.it)



DISTRIBUTED BY:



CE 2 YEAR WARRANTY